



Facebook Qrcode



Add us on Wechat




Scan here to shop



WhatsApp



ORIGIN DENTAL
SUPPLIES



ORIGIN DENTAL
SUPPLIES

ORIGIN DENTAL SUPPLIES

Head Office: 1/20 Barcoo Street, Roseville NSW 2069

Warehouse: 27/28 Barcoo St, Roseville NSW 2069

Phone: 02 9419 7658

Email: supply@origindental.com.au

Web: <https://www.origindental.com.au>

2023

FIRST EDITION

WWW.ORIGINDENTAL.COM.AU

Company Story

Origin Dental is a dentist-owned dental supply company. Setting up such a company is to leverage our existing buying power to import even a large quantity of low-cost consumables so that other smaller practices may also benefit from it. Origin Dental also offers a wide variety of other everyday consumables at a competitive price. All our products are tested by the dentists and service the dentists.

A COMPANY OF DENTISTS
TESTED BY DENTISTS
SERVICE DENTISTS

Origin Dental Supplies

Head Office: 1/20 Barcoo Street, Roseville NSW 2069
Warehouse: 27/28 Barcoo St, Roseville NSW 2069

<https://www.origindental.com.au>
02 9419 7658
supply@origindental.com.au



Facebook Qrcode



Add us on Wechat



Scan here to shop



WhatsApp

All pre orders will get 10% discount and items will be received in 20–30 days after each order

*T&C applies



Disposables	03–05
Endodontics	06–08
Infection Control	09–13
Impression	13–14
Restorative	15–19
Syringes and Needles	20
Phosphor Plates Readers	21–22
X– Ray	23
CBCT	24
3D Scanner	25–26
Handpieces	27–40
Bur	41–52
Trolley	53–54
Implant Motor	55–56
Oral Surgery	57–60
Hand Instrument	61–107
Othodontic	108–109
Other	110

Disinfecting Wipes Alcohol 18%

Cleans soil and contaminant from surfaces prior to disinfection. Medical/Dental grade. pH neutral, non-ionic. Non-staining and biodegradable. Non-corrosive to instruments and hand-pieces. Compatible with disinfection and sterilisation products. Approved for use with medical devices.



DP015
Packing: 220 pcs/bottle

Saliva Ejectors

Made of soft PVC, metal wire inside, easily forms into desired configuration, holds shape after bending, non-sterile.



DP012
Specification: clear tube with blue tip
Packing: 100 pcs/pack

Evacuation Tubes

Disposable, Softer end for patient comfort, no sharps, latex-free opaque plastic, high volume, non-sterile. Material : PE



DP013
Specification: white, non-vented "diagonal"
Packing: 100 pcs/bag

Disposable Air/Water Syringe Tips

Design with separated air and water channels, provides exceptional air/water spray.



DP001
Specification: Colorful, disposable
Packing: 250 pcs/bag

Dental Cups



DP004
Specification: White Paper Cups
Packing: 1000 pcs/carton

High Volume Suction Tips (Autoclavable)

Autoclavable High Volume Suction Tips (Blue)
No sharps, high volume, can be sterilised.



DP005
Specification: blue vented "S"
Packing: 50 pcs/bag

Disposable Mixing Wells

Best tools for adhesive, activators, and all material dispensed in small volume.



DP006
Specification: 2 Wells
Packing: 200 pcs/bag

Dental Bibs

Patient bibs are perfect for use during prophylaxis, cavity preparation and any other dental procedure where spray or splatter occurs. The tissue / poly design protects patient and equipment from moisture seeping through.



DP003
Specification: (13" x 18") (33x45cm)
2-ply paper (17g) + 1-ply PE (14.5g)
Color: Blue
Packing: 500 pcs/carton

Disposable Prophylaxis Brush

Premium quality prophylaxis brush, ideal for polishing of scaling operation, after filling and orthodontic.



DP010
Specification: Nylon
Packing: 100 pcs/box

Disposable Dappen Dishes

Fully disposable use, to save time and prevent cross-contamination for prophylaxis, etchant and bonding etc.



DP009
Specification: 23x32mm
Packing: 100 pcs/box

Aligners Chewies

Flexibility and elasticity
Reusable
No slippery
Suitable for all clear aligners and Invisalign



DP014
Specification: White, no flavor
Packing: 2 pcs/pack

Disposable Prophylaxis Cup

For low speed contra-angle handpiece with premium quality, and minimize chatter and heat.



DP011
Specification: Regular
Packing: 100 pcs/box

Cotton Rolls

100% cotton, for effective use in tight oral-cavity areas, non-linting, high absorbent, non-sterile.



DP002
Specification: #2, 1-1/2" × 3/8"
Packing: 2000 pcs/box

Surgical Tips



DP017
Packing: 20 pcs/box

Articulating paper 30um



DP018
Packing: 400 sheets/box

Non-Woven Gauze

Folded , low linting, excellent absorption capacity.
ply 30g / m²




DP007
Specification: 2" x2" (5x5cm)
Packing: 200 pcs/bag



DP008
Specification: 3" x3" (7.5x7.5cm)
Packing: 200 pcs/bag

N-Files (Rotay File)



- Offset mass of rotation.
- Unique saymmetrical rotary motion.
- Reduced engagement due to the swaggering effect which limits undesirable taper lock.
- Affords more cross-sectional space for enhanced cutting, loading,and augering debris.
- Greater flexibility than traditional NiTi sequence for navigating even challenging canals.
- Provides greater resistance to cyclic fatigue,reduce file separation risk.

Stainless Steel

Length(mm)	Size	Cross- Section	Packing
21/25/31	● X1、● X2、● X3		6pcs/packsingle or assortedsizes are available 3 pcs / pack if mixing size

TC-Files

- Imported CM-Wire NITi memory alloy Excellent fleibility & resistance to cyclic fatigue
- Triangular cross section deigu Keep center stiffness & high cutting efficiency
- Guiding non cutting force safety tip Efficient compliance with root canal guidance downward
- 04/06 Smaller tapered design Scientific sequence design



6 pcs / pack



0817 0315 0420 0425 0430

Stainless Steel


Model	Taper	Speed	Torque	Model	Taper	Speed	Torque
17#	.08	250rpm	3.0N.cn	25#	.04	250rpm	1.5N.cn
15#	.03	200rpm	0.8N.cn	30#	.04	300rpm	2.0N.cn
20#	.04	250rpm	1.0N.cn				

Open files



ENXA0817
Packing: 6 pcs/pack

K-Files



- For all preparation techniques with filling motion
- Search clean and cut dentine, make the wall smooth
- Square cross-section
- Stainless Steel/NITI
- Operation Type: Hand use/Engine use
- 6pcs per pack, available for single or assorted sizes

Stainless Steel/NITI

length	21mm、25mm、31mm
Model	06#、08#、10#、15#、20#、25#、30#

Endo Ruler (Blue)



740-601B
Packing: 1 pc/bag

Endo Irrigation Needle

Designed with a closed-end cannula and a side window for irrigation and safety in irrigating root canals and periodontal pockets.
EO sterile.




EN001
Specification: 27G, Grey
Packing: 100 pcs/box



EN014
Specification: 30G, Grey
Packing: 100 pcs/box

DTE Apex locator



EN015

An electronic apex locator is an electronic device used in endodontics to determine the position of the apical constriction and thus helps determine the precise working length.

- Using the latest multi-frequency apical positioning technology, Woodpecker DTEDPEXVApex Locator avoids the influence of subjective factors of measurement results.
- Measuring results are not affected by the tooth type, effects of root canal orifice size, or root canal internal environment.
- It shows a real-time display of the moving file in the root canal.

INDICATION:

- DTEDPEXVApex Locator is used in endodontics to determine working length precisely.

CONTRAINDICATION:

- The Woodpecker DTEDPEXVApex Locator is not recommended for use;
- In patients who have a pacemaker or other implanted electrical devices or have been cautioned by their physicians against the use of small electric appliances such as shavers, hairdryers, etc.
- In patients allergic to metals.
- In children.

ADVANTAGE:

- In clinical data matching, the accuracy rate has been improved to 97.71%.
- Upgraded chip, real-time response. The sensitivity has been improved to double.
- The innovative algorithm, multi-frequency measurement & anti-interference.
- LCD Screen, the file's trajectory can be clearly shown in different colours.
- Multi-frequency apical positioning technology.
- Avoiding the influence of subjective factors for measurement results.
- Measurement results will not be affected by tooth type.
- The Root canal's internal environment does not hamper its effectiveness.

Steps for using Woodpecker DTED pex V:

- When the file reaches the front region of the apical foramen, the white bars are displayed.
- When the file has gone to the position nearby the apical foramen, the green bars are displayed.
- The red bars are displayed when the file has exceeded the apical foramen. A continuous beep sound will be generated at the same time.

Barrier Film

Easily placed on any surface and removed without leaving any residue, reduces cross-contamination.




IC004
Packing: 1200 sheets/roll

Powder Free Latex Gloves



SKU	Specification	Packing
IC006	X-small	100 pcs/box, 10 boxes/carton
IC007	Small	100 pcs/box, 10 boxes/carton
IC008	Medium	100 pcs/box, 10 boxes/carton
IC009	Large	100 pcs/box, 10 boxes/carton

Plastic X-ray Digital Sensor Sleeve



IC003
Packing: 500 pcs/box

Autoclave Class 6 Test Strips

Placed within the autoclave with sterilisation pouches or instruments cassette, indicated whether the sterilant has penetrated the pack, changes colour when fully processed. Recommended frequency: every load.



IC001
Specification: class 6 test strips
Packing: 200 pcs/box

Powder Free Nitrile Gloves



SKU	Specification	Packing
IC010	Small	100 pcs/box, 10 boxes/carton
IC011	Medium	100 pcs/box, 10 boxes/carton
IC012	Large	100 pcs/box, 10 boxes/carton


Barrier Syringe Sleeve



IC002
Packing: 500 pcs/box

Self-Sealing Sterilization Pouches

Available for steam or E. O., with color changing indicator. All the sizes are the outer dimensions of the pouch.



SKU	Specification	Packing
IC014	Sterilization Self Sealing Pouch 135mm*260mm	200 pcs/bag
IC015	Sterilization Self Sealing Pouch 190mm*330mm	200 pcs/bag
IC016	Sterilization Self Sealing Pouch 57mm*100mm	200 pcs/bag
IC017	Sterilization Self Sealing Pouch 57mm*130mm	200 pcs/bag
IC018	Sterilization Self Sealing Pouch 70mm*230mm	200 pcs/bag
IC019	Sterilization Self Sealing Pouch 70mm*260mm	200 pcs/bag
IC021	Sterilization Self Sealing Pouch 90mm*165mm	200 pcs/bag
IC022	Sterilization Self Sealing Pouch 90mm*230mm	200 pcs/bag
IC023	Sterilization Self Sealing Pouch 90mm*260mm	200 pcs/bag
IC036	Sterilization Self Sealing Pouch 300mm*390mm	200 pcs/bag

Surgical Level 2 Face Mask

ARTG No : 340208
50PCS / BOX
Comfort-Fit
Splash Proof
No Latex
Fitting Compactness
Multi-layered Protection



IC013
Packing: 50 pcs/box

Infection Control

Bouffant cap

Bouffant cap (color:blue) , 21"(51cm), single elastic, 10g/square meter



IC033
Packing: 100 pcs/bag

Bibs Chain

Stretch tie with SS tip Colors mixing color



IC034
Packing: 1 pc/bag

Sterilization Tape Autoclave

Sterilization Tape Autoclave, 19mmX50m/rill



IC035
Packing: 1 pc/bag

Face Shield



IC027
Packing: 10 pcs/bag

Disposable CPE Gowns



IC025
Packing: 10 pcs/box

Disposable PP Gowns



IC026
Packing: 10 pcs/pack

Infection Control

Plastic Curing Light Tip Covers



IC031
Packing: 4.7"x1.2" 200 pcs/box

Sterilization Wrap



IC024
Packing: 60cm*60cm 500 pcs/carton

Plastic Curing Light Tip Covers



IC030
Packing: 4.33"x2.95" 200 pcs/box

Intraoral Camera Sheath



IC032
Packing: 500 pcs/bag

Headrest Cover

Disposable cover for headrest. Material:plastic.



IC028
Packing: 250 pcs/box
10"x14"
Thickness: 2s


Plastic Half Chair Covers

Disposable cover for backrest and headrest of chair.



IC029
Packing: 250 pcs/box
27-1/2"x24"
Thickness: 1.7s

Helix Indicators



IC037
5mm*73mm
Packing: 200 pcs/box

Impression

Mixing Tips

High performance , for 50ml dual-barrel impression material cartridges.



SKU	Specification	Packing
IM016	Yellow 5 . 0mm , (1 : 1) for impression material	50 pcs/bag
IM015	Green 6 . 5mm , (1 : 1) for impression material	50 pcs/bag
IM012	Blue / Orange , (10 : 1) for temporary crown material	50 pcs/bag
IM013	Brown / hite Short , (1 : 1) for cement	50 pcs/bag
IM014	For auto impression mixing systems , Dynamic Mixing Tips fits 3M PENRAMX II or II type mixing machine	50 pcs/bag
IM017	Kulzer Dynamix Mixing Tips 5:1 for Heraeus	50 pcs/bag

Dental Dispensing Gun

High performance impression gun fits all standard size cartridges.




IM010
Specification: 75ml, 1:1 / 2:1
Packing: 1 pc/box

IM011
Specification: 75ml, 4:1 / 10:1
Packing: 1 pc/box

Intra-Oral Tips

Tightly combination , easy dispensing.



IM009
Specification: Yellow, fit for IM-Mix-Yel-YMT
Packing: 100 pcs/bag

Impression Trays

Thick and wide edge trays , with high rigidity and toughness , not easy for demoulding, 100% latex-free.



SKU	Specification	Packing
IM007	Small Upper	12 pcs/bag
IM008	Small Lower	12 pcs/bag
IM005	Medium Upper	12 pcs/bag
IM006	Medium Lower	12 pcs/bag
IM003	Large Upper	12 pcs/bag
IM004	Large Lower	12 pcs/bag

3-In-1 Impression Trays



IM018
Specification: Sideless
Packing: 36 pcs/box

Retainer Boxes / Denture Baths

Sanitary type,
which also prevent to be spilled,
portable pocket box for denture.



IM002
Packing: 12 pcs/box

IM001
Packing: 12 pcs/box

Sectional Matrix System A2

The separating ring prevents the impaction.
The clamps with wedge-shaped tips separate teeth , fix positions and achieve tight contact;
The adaptive wedges and diamond wedges are contoured to help seal matrix at proximal box;
Dedicated instruments make the practice safe , easy and efficient.

Components can be purchased individually.



RE007

Specification: A2 Kit:- 3 rings / – 6 sizes of wedges / – 4 sizes of bands / – 3 sizes wedge guards / – 2 forceps per box

Clamp with SmallWedaē-shaped Tips



RE009



RE010



Bicuspids



Molar



Tall molar



Molar with extensions

Diamond Wedges



Extra small, small, medium, large, deep

Adaptive Wedges



Standard

Wedging Guards



small, medium, large

Forceps for placing sectional matrix bands



791-223

Forceps for placing rings



791-224


Adaptive Wedges

New design , Inner part is hard plastic , outer part is soft silicone , which make the balance between hardness and elastic.
The elliptical inner core and silicone can be match with various adjacent interdental spaces , and press the matrix bands on the teeth closely.
Autoclavable.



RE021

Specification: Plastic core, covered with silicone
Packing: 40 pcs/box



Diamond Wedges

New design Ideal contour to seal matrix / tooth interface.
Collapsing tip for easy placement , it rebounds when fully seated.
Autoclavable.



SKU	Specification	Packing
RE017	Extra Small	40 pcs/box
RE020	Small	40 pcs/box
RE016	Medium	40 pcs/box
RE015	Large	40 pcs/box
RE019	Deep	40 pcs/box
RE018	Longer Medium	40 pcs/box
RE014	Large Green	40 pcs/box

Matrix Bands

Tight contacts , excellent retention & separating teeth without impaction.



RE022
Bicuspid




RE023
Molar



RE024
Tall molar



RE025
Molar with extensions



RE026
4 Size Sectional Matrix Bands (15/each)
Packing: 60 pcs/box

RE027
4 Size Sectional Matrix Bands (25/each)
Packing: 100 pcs/box


Shading Dics



RE006
Packing: 1 pcs/bag


Micro Applicators

Ideal for cavity linings , bonding agents , sealants , etchants , hemostatic solutions.
Non-absorbent , bendable neck , non-linting.
Material : PP



SKU	Specification	Packing
RE005	Supper-fine, White	100 pcs/pack, 4 packs/box
RE004	Fine, Mauve	100 pcs/pack, 4 packs/box
RE003	Regular, Blue	100 pcs/pack, 4 packs/box

Mixing pad



RE029
Packing: 7.6*7.6cm 50 pcs/book

RE030
Packing: 5.1*5.1cm 50 pcs/book

O-light II


More-focused light output

Wide spectrum
high light intensity


RE033
Packing: 1 pc/bag

Upgraded focused light output

Unique Woodpecker technology of focused light output reduces 55.87% light angle, allowing more concentrated energy and higher curing efficiency.



O-Light II with US A brand



O-Light II with US B brand



Magnetic Point Cure Lens

Provides pinpoint curing of small composites and is helpful for tack curing veneers and all porcelain crowns.

3 modes for option



Normal mode
Intensity
1000-1200 mw/cm²
Time setting
5s, 10s, 15s, 20s



Normal mode
Intensity
1800-2000 mw/cm²
Time setting
3s, 5s




Normal mode
Intensity
2700-3000 mw/cm²
Time setting
1s, 3s

Upgraded wide spectrum

The emission wavelength range of O-light II Wide Spectrum Curing Light is 385nm-515nm, mainly adding the waveband of 385nm-420nm. The additional band can better curing the resin materials containing new-style photoinitiators such as TPO. For example, restoration materials such as Tetric N-ceram Bulk Fill (Ivoclar), Amelogen Plus (Ultradent), and ENA Hn (Micerium). The additional wave band can also reduce the risk of resin shedding and incomplete curing.

O-light II Spectral Curve and Resin Spectral Absorptivity



The light absorbing peak of traditional materials (containing camphorquinone, the red line above) is 460nm, so the curing can be done both ordinary Curing Light and Wide Spectrum Curing Light. However, the camphorquinone is yellowish and it is not suitable to add a large amount of it to bleached resin or light-colored resin. The light absorbing peak of new material (including TPO, the blue line above) is 405nm. O-light II owns an additional wave band of 385nm-420nm which can better cure these materials. TPO is white and is suitable for being added to bleached resins and light-colored resin.

1 SEC FOR CURING

Capsule Activator GUN



RE002
Packing: 1 pcs/box

Composite Gun Metal Tip


Ergonomically designed dispensing gun used with most unidos tips.



RE001
Packing: 1 pc/bag


Capsule Dispenser Gun (Autoclavable)

Ergonomically designed dispensing gun used with most unidos tips. Autoclavable up to 134°C.



RE008
Packing: 1 pc/bag

Tofflemire Matrix Bands



RE012
Packing: #1 Ultra Thin 0.0015 144 pcs/pack
RE013
Packing: #2 Ultra Thin 0.002 144 pcs/pack

Curing Light

Digital display
Curing adjustable
Membrane key



RE028
Packing: 1 pc/box

Monojet Irrigation Syringe

Disposable, non-sterile syringes ideals for endo irrigation or intra-oral impression delivery.



SN001
Specification: Clear
Packing: 50 pcs/bag

Dispenser Tips

Disposable blunt end pre-bent needle tips for use with etchants, resins and flowable composites, non-sterile.




SN002
Specification: 19G
Packing: 100 pcs/bag



SN003
Specification: 22G
Packing: 100 pcs/bag

Dental Needle



SN004
Packing: 30G*25mm 100 pcs/bag

Phosphor Plates Readers

DTE Free Scan



- Ultra-high image resolution provides physicians with clearer diagnostic images
- High-quality user interface makes shooting and reading easier
- High fidelity image storage and transmission, easy to refer to and follow-up diagnosis
- Fast reliable scanning processing provides the excellent user experience



Phosphor Plates Readers

Easy to operate


One-button start for scanning, 120-second auto lock to keep datasafe.

Connecting to a computer to achieve direct operation without the need for dental software. This method can help simplify the workflow and improve efficiency.



Advanced algorithm leads to intelligence

Advanced algorithm leads to intelligence
Multi-layered structural image algorithm framework enhances the details of super-resolution images
Multi-mode processing algorithm to meet the image viewing habits of different doctors



Quick image viewing

Imported ultra-thin imaging plate has moderate thickness to keep oral cavity comfortable.
The 0-2 size imaging plates are available for photographing all tooth positions.



0# XR004 1# XR005 2# XR006
Packing: 1 pc/bag

Powerful chip Intelligent scanning

FPGA chip accelerates data acquisition.
Support image viewing on mobilephones, tablets, and computers.
Exclusive image viewing software.
Greatly facilitate image viewing and processing.



Real-time intelligent dose analysis Safe dental film photographing

Appropriate dose.
Accurately record the dose value for each photographing.
Timely reminder of exposure dose for clear images.
Ensure secure photographing.



Good-looking and handy Compacted size of 0.046m² and space-saving position

The silver-black scheme is suitable for upscale and official occasions
The compacted size of 0.046m² remains full-fledged
Standing nearby armchair to safeguard your teeth

Barrier Envelopes

Strong seal protects the integrity of the phosphor plate, with easy teared PE film, soft, easy loading and opening / removing.
Mouth way on two neck parts for easy tearing.
Widen design on two neck parts ensures the edge not be penetrated when doing edge oblique mouth way , better protects the digital phosphor plate from cross-contamination.



XR001
Model: SIZE 0
Size: 22 × 25mm
Packing: 300 pcs/box



XR002
Model: SIZE 1
Size: 24 × 40mm
Packing: 300 pcs/box



XR003
Model: SIZE 2
Size: 31 × 41mm
Packing: 300 pcs/box

Medical Imaging

Product Introduction
Dental Imaging Equipment—cbct

The Meyer extraoral CBCT equipment collects full oral data in one scan and reconstructs all aspects of high-resolution images as needed for accurate clinical diagnostics . The resulting 3D images and analytical data provide essential basis for dental filling , implant and orthodontics

The product name		Meyer CBCT—Wise	Meyer CBCT—Honor	Meyer CBCT—Dream
Maximum FOV (cm)		12×10	15×10.5	17×11
Clinical Advantages		1.Full scan from lower jaw to maxillary sinus. 2.Suitable for general + local diagnostics, and dental implant.	1.Full scan from lower jaw to maxillary sinus and TMJs. 2.Suitable for general + local diagnostics, TMJs examination and dental implant.	1.Full scan from lower jaw to maxillary sinus and TMJs, high resolution part CT. 2.Suitable for general + local diagnostics, TMJs examination and dental implant.
Technical Parameters	Tube Voltage	60~90kV		
	Tube Current	2~10mA		
	X ray radiation time	CT~20secs Pano~17secs Ceph-lateral~12secs Local CT~11.75secs		
	Focal spot to skin	250mm		
	Power supply	AC 220V 16A		
	Tube Model	D-054SB		
	Dimension	3DE High: 1755mm, Wide:2485mm, Length:1415mm 3D High: 1755mm, Wide:1090mm, Length:1415mm		
	Net Weight	3DE: 260KG 3D: 220KG		
Operating System		Windows 7 or above, Intel Core i5 or above, 32-bit/64-bit		



SS-X9010DPRO-3D / 3DE

Aoralscan3 Intra Oral Scanner



Scan Field	Standardscannertip: 16mm × 12mm × 22mm Miniscannertip: 12mm × 9mm × 22mm
Scan Depth	–2–20mm from exit surface of tip
Scan Principle	Non-contact scanner with structured light
Dimension (L x W x H)	Dimension(L × W × H)285mm × 33mmx46mm
Weight	240 ± 10g(without cables)
Output	STL, OBJ, PLY
Connector	USB 3.0
Power Input	12VDC/3A


Recommended PC Configuration

CPU	Intel Core i 7-8700or higher
Memory	16GBor higher
HardDisk Drive	256GBSSDor above
Graphic Card(GPU)	NVIDIA® RTX20606GBor higher
OperatingSystem	Windows 10 Professional(64-bit) or later versions of Windows operatingsystem
Display Resolution	1920 × 1080, 60Hzor higher
I/O Ports	More than 2 type-A USB 3.0(or higher) ports

Empowering a Highly Efficient Digital Chairside Solution

Say hello to SHINING 3D 's brandnewintraoral scanner, A oral scan 3!Exceptional, smart and elegant, A oral scan 3 improves the overall clinical experience.It allows users to easily acquire vivid digital impressions in a faster and more accurate way.

Moreover, A oral scan 3 comes with a wide range of different intelligent functions, guar-an teeing more comfortable chairside experi-ence and efficient clinic-lab collaboration.



User-Friendly
Smart and Intuitive–Smooth Scanning Experience for Dentists


With powerful and intelligent data processing cap abl ities, Aoralscan3is easy to maneuver.Even begin-ners can get ideals can results in minutes.



New UI
Neater and more user-friendly interface
More flexible workflow



AI Scan
Automatically identifies and filters out unnecessary soft tissue data during scanning resulting in a quicker and cleaners can process.



Motion Sensing
Dentists can scan with minimal contact with the computer, ensuring a safe and efficient scan experience.



Cart Option
The dedicated cart solution was designed to professionally round off the Aoralscan experience in a broad array of working scenarios.

Dental Air Polishing Handpiece

AP-H Plus

Comfortable Clean Lightweight
Easy Cleaning Mode

HP-AP01

replaceable tip
HP-AP02



Supragingival & subgingival mode
Two-in-one



0.7mm small caliber nozzle



Three-section air-polishing
handpiece design



360° rotatable
Pre-bendable
subgingival nozzle



Lateral transparent powder
container enables easier
observation of remaining powder



The whole device
weights about 125g




HP 1:3

Woodpecker 45° Hand Dental 1:3 Contra Angle LED Handpiece External Water

Speed ratio of 1:3
Autoclavable design
Improved design for extra comfort
German Bearing for superior life and performance
Model : WJ-SG45CL
One year warranty

Handpieces



HP-HW5L

Dental LED Light Ultrasonic Scaler Handpiece
HW-5L For UDS LED Scaler

HW-5L LED Handpiece:
Dental ultrasonic scaler handpiece
Autoclavable under the high temperature of 135°C and pressure of 0.22Mpa
With 6 LED
compatible with EMS
fit for Woodpecker UDS series ultrasonic scalers like
UDS-E LED, UDS-N2 LED, UDS-N3 LED, UDS-K LED, UDS-P LED

Handpieces

J4 Mini head dental handpiece for children

Mini head high speed handpiece, NSK type quick coupling, 4 holes, equipped with short shank burs.



J4-MUQ-M4

Buy One get One FREE! FREE Coupling!

8.9mm Extraordinary miniature head

Every dentist deserves to have a piece of mini head dental handpiece for children -- J4



J4 Mini head dental handpiece for children



Standard dental handpiece

Specifications

Model	MUQ
Product code	J4
Order number	100068
Working pressure	0.23MPa-0.25MPa
Rotating speed	380,000rpm-450,000rpm
Noise	≤60dB
Atomized water pressure	0.3-0.35MPa
Working water pressure	0.2-0.23MPa
Water flow	90-110ml/min
Radial vibration	Vibration value for standard testing bar≤0.03mm
Surface	Titanium plating
Water spray	single spot
Chucking type	Push-button type
Cartidge	Exposed cartridge
Bearing	German ceramic bearing
Autoclave	Can be sterilized no more than 135 degree
Connection	4 holes quick coupling, compatible with NSK

J6K&J6N Fiber Optic Handpiece

Quiet noise , Big torque , Long lite span , Strong illumination



J6K



K Coupling


J6 Fiber Optic Handpiece

Quiet noise , Big torque , Long lite span , Strong illumination.



Consistent electric current flow

There is no impact to illumination of light no matter what voltage is used in the dental chairs. it can even bright up to maxLUX at IV.



Polarity-freeIt

It can be used as long as it connects to any brand of dental chairs without paying attention to whether positive and negative pole is consistent with those in dental units.



Only 46g

Only 46g in weigh

Decrease feeling of fatigue during long time grip for use.



German bearings

German bearings, 1 year warranty.



FREE Coupling!

Fiber optic high speed handpiece, Nsk type quick coupling, 6 holes.

J6-TUQL-N-M6

Special Deal
**Buy One get
One FREE!**

Handpieces



Which is more practical compared
to air-driven handpieces

Digital rotation speed could have a real-time view to make the treatment more precise

Handpieces

Special Deal
Buy One get One FREE!

1:1 internal water spray contra angle with fiber optic



1:1 L

Rotation ratio	1: 1
Connector	standard E-type connector
Rotating speed	18000rpm-20000rpm
Noise	≤60dB
Light	Build-in LED
Water spray	single spot
Chucking type	Push-button type
Product model	101113

1:5 electric contra angle with fiber optic



1:5 L

Rotation ratio	1: 5
Connector	standard E-type connector
Rotating speed	highest 200000rpm
Noise	≤60dB
Light	triple spray
Water spray	optical fiber lighting
Chucking type	Push-button type
Product model	105019

Handpieces



It is more convenient to operate with internal water spray instead of water pipe claws connected.

Handpieces

M2 Internal water spray low speed handpiece

Internal Watercourse Straight Handpiece.

Special Deal
Buy One get One FREE!



M2-SH



Push button type



Fixed spray hole



Water spray



Internal watercourse

Specifications

Model	LN
Product code	M2
Order number	102013
Working pressure	0.28MPa-0.3MPa
Rotating speed	1800-20000rpm
Noise	≤65dB
Atomized water pressure	0.3-0.35MPa
Working water pressure	0.2-0.23MPa
Water flow	90-110ml/min
Surface	chrome plating
Water spray	Internal water spray
Chucking type	push button type for contra angle, tension ring for straight head
Bearing	ceramic bearing
Autoclave	Can be sterilized no more than 135 degree
Connection	2 holes (Borden)/4 holes (Midwest) coupling



J45

LED minimally invasive dental handpiece with light source

Essential tool for various teeth minimally invasive operations

J45 Degree Surgical Handpiece

45 degree, easily for the back teeth surgical, provides a very good angle to directly deliver air, water or mist to the mouth, which is convenient for operation, doctors may not feel tired after long operation.

Special Deal
Buy One get One FREE!



J45-LED



100% rear-exhaustion



With LED



With special tungsten steel bur



limited vibration

Specifications

Model	TU
Product code	J45
Order number	100070
Working pressure	0.25MPa-0.3MPa
Rotating speed	360,000rpm-420,000rpm
Noise	≤60dB
Atomized water pressure	0.3-0.35MPa
Working water pressure	0.2-0.23MPa
Water flow	90-110ml/min
Radial vibration	Vibration value for standard testing bar ≤0.03mm
Surface	Titanium plating
Water spray	single spot
Chucking type	Push-button type
Cartidge	Exposed cartridge
Bearing	German ceramic bearing
Autoclave	Can be sterilized no more than 135 degree
Connection	optional LED

Handpieces



Easy to disassemble , high torque , long life span

Handpieces

20 : 1 implant contra-angle handpiece



Easy to disassemble
The perfect structure of easily assembling and disassembling is conducive to thoroughly clean and sterilize.



Little vibration
The German bearings can rotate very smoothly, its radial runout which is less than 0.03mm, is equal to a human hair, meeting all kinds of precision Operation.



High torque
Due to the unique structure of gear, its max torque can be 80N cm, perfectly meet the needs of dental implant treatment as much as possible.

Common accessories



Contra angle middle gear
NO : 210128



Cartridge
NO : 210126



Back cap
NO : 200182



Blue oil seal
NO : 316025


Technical Parameters

Max T orgue	80N · cm
Noise	≤50dB
Radial runout	≤0.03mm
Illumination	LED Fiber Optic
Gear Ratio	20 : 1 Reductiono.
Head diameter	Φ9.6mm
Chucking type	push button
Sterilization	No more than 135℃
Order number	112032

Diamond Bur

Advantages

- 1. Strictly comply with ISO Standard
- 2. Sharper and More Efficient
- 3. Durable and Valuable
- 4. Full Range of shapes and grits
- 5. Fastly develop new items
- 6. Individual Design and OEM Service








- Grit
- black / Super Coarse
 - Green / Coarse
 - Blue / Medium
 - Red / Fine
 - Yellow / Extra Fine

Standard




 0.8 19.0 001(801) BR-49	 1.0 19.0 001(801) BR-45	 1.2 19.0 001(801) BR-46	 1.4 19.0 001(801) BR-41	 1.6 19.0 001(801) BR-40	 1.8 19.0 001(801) BR-31	 2.1 19.0 001(801) BR-28	 2.3 19.0 001(801) BR-26	 2.5 19.0 001(801) BR-25	 2.9 19.0 001(801) BR-23
001-008M 001-010M 001-012M 001-014M 001-016M 001-018M 001-021M 001-023M 001-025M 001-029M									
Ball, Round / BR									
 * 1.2 22.0 697(801L) BR-L46	 * 1.4 22.0 697(801L) BR-L41	 * 1.6 22.0 697(801L) BR-L40	 * 1.8 22.0 697(801L) BR-L31	 1.2/2.3 19.0 002(802) BC-43	 1.4/2.3 19.0 002(802) BC-42	 1.6/3.5 19.0 002(802) BC-32	 1.8/2.3 19.0 002(802) BC-31	 1.0(0.8)/1.0 19.0 010(805) SI-45	 1.2(0.9)/1.2 19.0 010(805) SI-46
697-012M 697-014M 697-016M 697-018M 002-012M 002-014M 002-016M 002-018M 010-010M 010-012M									
Round with Long Neck / BR Round with Collor / BC Inverted Conical / SI									
 1.4(1.0)/1.6 19.0 010(805) SI-47	 1.6(1.1)/1.8 19.0 010(805) SI-48	 1.0/2.3 19.0 019(806) DI-41	 1.4/3.2 19.0 019(806) DI-42	 3.2/5.0 19.0 039(811L) EX-11	 3.5/7.0 20.0 039(811L) EX-12	 4.5/0.9 19.0 041(818) WR-44	 4.2/1.8 19.0 068(909) WR-13	 1.0/4.0 19.0 109(836) SF-41	 1.2/4.0 19.0 109(836) SF-31
010-014M 010-016M 019-010M 019-014M 039-032M 039-035M 041-045M 068-042M 109-010M 109-012M									
Double Conocal Wheel Flat End Cylinder									
 1.4/4.0 19.0 109(836) SF-42	 1.4/6.0 21.0 110(837) SF-21	 1.6/6.0 21.0 110(837) SF-15	 1.2/8.0 22.0 111(837L) SF-11	 1.4/8.0 22.0 111(837L) SF-12	 1.6/8.0 22.0 111(837L) SF-13	 1.2/8.0 22.0 130(885) CP-11	 1.4/8.0 22.0 130(885) CP-12	 1.0/4.0 19.0 139(838) SR-41	 1.2/4.0 19.0 139(838) SF-31
109-014M 110-014M 110-016M 111-012M 111-014M 111-016M 130-012M 130-014M 139-010M 139-012M									
Cylinder, end cone, pointed Round End, Cylinder									
 1.4/4.0 19.0 139(838) SR-42	 1.4/6.0 21.0 140(880) SR-21	 1.2/8.0 22.0 141(881) SR-11	 1.4/8.0 22.0 141(881) SR-12	 1.6/8.0 22.0 141(881) SR-13	 1.2/0.2 21.0 150(839) EX-18	 1.2/6.0 21.0 164(852) TC-26	 1.4/8.0 22.0 165(858) TC-21	 1.0/10.0 23.0 166(859) TC-27	 1.2/10.0 23.0 166(859) TC-28
139-014M 140-014M 141-012M 141-014M 141-016M 150-012M 164-012M 165-014M 166-010M 166-012M									
End Cut					Needle				

 1.4/10.0 23.0 166(859) TC-10	 1.6/10.0 23.0 166(859) TC-11	 1.8/10.0 23.0 166(859) TC-16	 1.0(0.6)/3.5 19.0 169(845) TF-41	 1.2(0.7)/4.0 19.0 170(845) TF-42	 1.4(0.8)/4.0 19.0 170(845) TF-43	 1.6(0.9)/4.0 19.0 170(845) TF-31	 1.8(1.0)/5.0 20.0 170(845) TF-23	 2.1(1.2)/5.0 20.0 170(845) TF-22	 1.4(0.7)/6.5 21.0 171(846) TF-20
166-014M	166-016M	166-018M	169-010M	170-012M	170-014M	170-016M	170-018M	170-021M	171-014M
Flat End Cone									
 1.6(0.8)/6.5 21.0 171(846) TF-21	 1.0(0.6)/8.0 22.0 172(847) TF-07	 1.2(0.7)/8.0 22.0 172(847) TF-08	 1.4(0.7)/8.0 22.0 172(847) TF-09	 1.6(0.8)/8.0 22.0 172(847) TF-10	 1.8(0.9)/8.0 22.0 172(847) TF-56	 2.3(1.3)/9.0 22.0 172(847) TF-14	 1.4(0.7)/10.0 24.0 173(848) TF-11	 1.6(0.8)/10.0 22.5 173(848) TF-12	 1.8(0.9)/10.0 22.5 173(848) TF-13
171-016M	172-010M	172-012M	172-014M	172-016M	172-018M	172-023M	173-014M	173-016M	173-018M
 1.0(0.6)/4.0 19.0 196(855) CR-41	 1.2(0.7)/4.0 19.0 196(855) CR-31	 1.4(0.8)/4.0 19.0 196(855) CR-42	 1.6(0.9)/4.0 19.0 196(855) CR-12	 2.0(1.0)/5.0 20.0 196(855) CR-11	 1.4(0.7)/6.5 21.0 197(849) TR-20	 1.6(0.8)/6.5 21.0 197(849) TR-21	 1.8(1.0)/6.5 21.0 197(849) TR-62	 2.0(1.1)/6.5 21.0 197(849) TR-63	 1.4(0.8)/8.0 22.0 198(850) TR-25
196-010M	196-012M	196-014M	196-016M	196-020M	197-014M	197-016M	197-018M	197-020M	198-014M
Round End Cone									
 1.6(0.9)/8.0 22.0 198(850) TR-24	 1.8(1.0)/8.0 22.0 198(850) TR-23	 1.4(0.7)/10.0 22.0 199(852) TR-11	 1.6(0.8)/10.0 22.0 199(852) TR-12	 1.8(0.9)/9.0 22.0 199(852) TR-13	 2.2(1.3)/9.0 22.0 199(852) TR-14	 2.2(1.3)/11.0 24.0 200(852L) TR-15	 2.4(1.4)/11.0 24.0 200(852L) TR-19	 1.0(0.8)/2.8 19.0 233(808) EX-41	 1.4(1.1)/4.0 19.0 234(808L) EX-31
198-016M	198-018M	199-014M	199-016M	199-018M	199-022M	200-022M	200-024M	233-010M	234-014M
Pear									
 1.8(1.2)/5.0 19.0 237(808R) EX-20	 2.1(1.6)/5.0 19.0 237(808R) EX-21	 3.2(1.9)/5.0 19.0 237(808R) EX-26	 1.2/6.0 21.0 248(861) FO-41	 1.4/6.0 21.0 248(861) FO-20	 1.4/8.0 22.0 249(862) FO-21	 1.6/8.0 22.0 249(862) FO-22	 1.2/10.0 24.0 250(863) FO-11	 1.4/10.0 24.0 250(863) FO-12	 1.6/10.0 24.0 250(863) FO-13
237-018M	237-021M	237-032M	248-012M	248-014M	249-014M	249-016M	250-012M	250-014M	250-016M
Flame									
 1.8/3.6 19.0 257(368,830) FO-32	 2.3/5.0 19.0 257(368,830) FO-34	 2.8/5.0 19.0 257(368,830) FO-25	 3.2/5.0 19.0 257(368,830) FO-27	 2.0/4.0 19.0 277(379,833) TR-22	 3.0/5.0 19.0 277(379,833) TR-28	 1.2/6.0 21.0 288(867,877) SO-20	 1.4/8.0 22.0 289(868,878) SO-21	 1.6(1.0)/8.0 22.0 298(878K) 878K16	 1.8(1.1)/8.0 22.0 298(878K) 878K18
257-018M	257-023M	257-028M	257-032M	277-020M	277-030M	288-012M	289-014M	298-016M	298-018M
Egg									
Torpedo									
 1.6(1.0)/10.0 24.0 299(879K) 879K16	 1.4/8.0 22.0 115(842SP) BC-11	 1.6/8.0 22.0 178(848SP) BC-12	 cutting depth 0.3mm 1.6(1.0)/6.0 21.0 552(834) VP-20	 cutting depth 0.4mm 1.8(1.0)/6.0 21.0 552(834) VP-21	 cutting depth 0.5mm 2.1(1.1)/6.0 21.0 552(834) VP-22	 1.6(0.8)/8.0 22.0 546(847KR) 847KR16	 1.8(0.9)/10.0 24.0 553(848KR) 848KR18		
299-016M	115-014M	178-016M	552-016M	552-018M	552-021M	546-016M	553-018M		
Depth Marking								Taper, Round Edge	

Super Coarse Black

 Super Coarse Black	 1.4 19.0 001(801) BR-41S	 1.6 19.0 001(801) BR-40S	 1.8 19.0 001(801) BR-31S	 2.1 19.0 001(801) BR-28S	 2.3 19.0 001(801) BR-26S	 2.5 19.0 001(801) BR-25S	 2.9 19.0 001(801) BR-23S	 1.6(1.1)/1.8 19.0 010(805) SI-48S	 3.2/5.0 19.0 039(811L) EX-11S
	001-014S	001-016S	001-018S	001-021S	001-023S	001-025S	001-029S	010-016S	039-032S
 3.5/7.0 20.0 039(811L) EX-12S	 1.2/4.0 19.0 109(836) SF-31S	 1.4/4.0 19.0 109(836) SF-42S	 1.2/8.0 22.0 111(837L) SF-11S	 1.4/8.0 22.0 111(837L) SF-12S	 1.6/8.0 22.0 111(837L) SF-13S	 1.6/10.0 23.0 166(859) TC-11S	 1.8/10.0 23.0 166(859) TC-16S	 1.6(0.8)/6.5 21.0 171(846) TF-21S	 1.4(0.7)/10.0 24.0 173(848) TF-11S
039-035S	109-012S	109-014S	111-012S	111-014S	111-016S	166-016S	166-018S	171-016S	173-014S
 1.6(0.8)/10.0 22.5 173(848) TF-12S	 1.8(0.9)/10.0 22.5 173(848) TF-13S	 1.6(0.8)/6.5 21.0 197(849) TR-21S	 1.4(0.7)/10.0 22.0 199(852) TR-11S	 1.6(0.8)/10.0 22.0 199(852) TR-12S	 1.8(0.9)/9.0 22.0 199(852) TR-13S	 1.8(1.2)/5.0 19.0 237(808R) EX-20S	 2.0/4.0 19.0 277(379,833) TR-22S		
173-016S	173-018S	197-016S	199-014S	199-016S	199-018S	237-018S	277-020S		

Coarse Green





 Coarse Green	 1.4 19.0 001(801) BR-41C	 1.6 19.0 001(801) BR-40C	 1.8 19.0 001(801) BR-31C	 2.1 19.0 001(801) BR-28C	 2.3 19.0 001(801) BR-26C	 2.5 19.0 001(801) BR-25C	 2.9 19.0 001(801) BR-23C	 1.6(1.1)/1.8 19.0 010(805) SI-48C	 4.2/1.8 19.0 068(909) WR-13C
	001-014C	001-016C	001-018C	001-021C	001-023C	001-025C	001-029C	010-016C	068-042C
 1.4/6.0 21.0 110(837) SF-21C	 1.4/8.0 22.0 165(858) TC-21C	 1.8(0.9)/10.0 22.5 173(848) TF-13C	 1.8/9.5 21.0 190 FO-54C	 1.8(1.0)/6.5 21.0 197(849) TR-62C	 1.8(0.9)/9.0 22.0 199(852) TR-13C	 1.8/3.6 19.0 257(368,830) FO-32C			
110-014C	165-014C	173-018C	190-018C	197-018C	199-018C	257-018C			

Fine Red

 Fine Red	 1.2 19.0 001(801) BR-46F	 1.4 19.0 001(801) BR-41F	 1.6 19.0 001(801) BR-40F	 1.8 19.0 001(801) BR-31F	 2.1 19.0 001(801) BR-28F	 2.3 19.0 001(801) BR-26F	 1.4/6.0 21.0 110(837) SF-21F	 1.2/0.2 21.0 150(839) EX-18F	 1.4/8.0 22.0 165(858) TC-21F
	001-012F	001-014F	001-016F	001-018F	001-021F	001-023F	110-014F	150-012F	165-014F

 1.6/10.0 23.0 166(859) TC-11F	 1.6(0.8)/6.5 21.0 171(846) TF-21F	 1.4(0.7)/10.0 24.0 173(848) TF-11F	 1.6(0.8)/10.0 22.5 173(848) TF-12F	 1.8(0.9)/10.0 22.5 173(848) TF-13F	 1.6(0.9)/4.0 19.0 196(855) CR-12F	 2.0(1.0)/5.0 20.0 196(855) CR-11F	 1.4(0.7)/6.5 21.0 197(849) TR-20F	 1.6(0.8)/6.5 21.0 197(849) TR-21F	 1.8(1.0)/6.5 21.0 197(849) TR-62F
166-016F	171-016F	173-014F	173-016F	173-018F	196-016F	196-020F	197-014F	197-016F	197-018F
 1.4(0.7)/10.0 22.0 199(852) TR-11F	 1.6(0.8)/10.0 22.0 199(852) TR-12F	 1.8(0.9)/9.0 22.0 199(852) TR-13F	 1.4(1.1)/4.0 19.0 234(808L) EX-31F	 2.1(1.6)/5.0 19.0 237(808R) EX-21F	 1.4/8.0 22.0 249(862) FO-21F	 1.6/8.0 22.0 249(862) FO-22F	 1.6/3.2 19.0 254(831) FO-30F	 1.4/8.0 22.0 289(868,878) SO-21F	
199-014F	199-016F	199-018F	234-014F	237-021F	249-014F	249-016F	254-016F	289-014F	

Extra Fina Yellow

 Extra Fina Yellow	 1.2 19.0 001(801) BR-46EF	 1.4 19.0 001(801) BR-41EF	 1.6 19.0 001(801) BR-40EF	 1.8 19.0 001(801) BR-31EF	 2.1 19.0 001(801) BR-28EF	 2.3 19.0 001(801) BR-26EF	 1.4/8.0 22.0 165(858) TC-21EF	 1.0/10.0 23.0 166(859) TC-27EF	 1.2/10.0 23.0 166(859) TC-28EF
	001-012EF	001-014EF	001-016EF	001-018EF	001-021EF	001-023EF	165-014EF	166-010EF	166-012EF
 1.6/10.0 23.0 166(859) TC-11EF	 1.6(0.8)/6.5 21.0 171(846) TF-21EF	 1.4(0.7)/10.0 24.0 173(848) TF-11EF	 1.6(0.8)/10.0 22.5 173(848) TF-12EF	 1.8(0.9)/10.0 22.5 173(848) TF-13EF	 1.6(0.9)/4.0 19.0 196(855) CR-12EF	 2.0(1.0)/5.0 20.0 196(855) CR-11EF	 1.6(0.8)/6.5 21.0 197(849) TR-21EF	 1.8(1.0)/6.5 21.0 197(849) TR-62EF	 1.6(0.9)/8.0 22.0 198(850) TR-24EF
166-016EF	171-016EF	173-014EF	173-016EF	173-018EF	196-016EF	196-020EF	197-016EF	197-018EF	198-016EF
 1.4(0.7)/10.0 22.0 199(852) TR-11EF	 1.6(0.8)/10.0 22.0 199(852) TR-12EF	 1.8(0.9)/9.0 22.0 199(852) TR-13EF	 2.1(1.6)/5.0 19.0 237(808R) EX-21EF	 1.0/4.3 23.5 243(831) FO-40EF	 1.2/6.0 21.0 248(861) FO-41EF	 1.4/6.0 21.0 248(861) FO-20EF	 1.4/8.0 22.0 249(862) FO-21EF	 1.6/8.0 22.0 249(862) FO-22EF	 1.6/3.2 19.0 254(831) FO-30EF
199-014EF	199-016EF	199-018EF	237-021EF	243-010EF	248-012EF	248-014EF	249-014EF	249-016EF	254-016EF
 2.3/5.0 19.0 257(368,830) FO-34EF									
257-023EF									

GF 313 Short 16.5mm

 1.0 16.5 001(801) BR-S45	 1.4 16.5 001(801) BR-S41	 1.2(0.9)/1.2 16.5 010(805) SI-S46	 1.6(1.1)/1.8 16.5 010(805) SI-S48	 1.4/8.0 19.0 165(858) TC-S21	 1.0(0.6)/3.5 16.5 169(845) TF-S41	 1.6(0.9)/4.0 16.5 170(845) TF-S31	 1.8(1.0)/5.0 17.5 170(845) TF-S23	 1.6(0.8)/6.5 18.0 171(846) TF-S21	 1.4(0.7)/6.5 18.0 171(846) TF-S20
001-S010M	001-S014M	010-S012M	010-S016M	165-S014M	169-S010M	170-S016M	170-S018M	171-S016M	171-S014M
 1.6(0.8)/6.5 18.0 197(849) TR-S21	 1.4/8.0 19.0 249(862) FO-S21								
197-S016M	249-S014M								

FG 316 Surgical 25mm

 1.0 25.0 001(801) BR-X45	 1.2 25.0 001(801) BR-X46	 1.4 25.0 001(801) BR-X41	 1.6 25.0 001(801) BR-X40	 1.8 25.0 001(801) BR-X31	 1.2(0.9)/1.2 25.0 010(805) SI-X46	 1.2/8.0 28.0 111(837L) SF-X11	 1.4/8.0 28.0 111(837L) SF-X12	 1.6/8.0 28.0 111(837L) SF-X13	 1.6/10.0 29.0 166(859) TC-X11
001-X010M	001-X012M	001-X014M	001-X016M	001-X018M	010-X012M	111-X012M	111-X014M	111-X016M	166-X016M
 1.4(0.7)/10.0 29.0 173(848) TF-X11	 1.6(0.8)/10.0 28.5 173(848) TF-X12	 1.8(0.9)/10.0 28.5 173(848) TF-X13	 1.6(0.8)/10.0 28.0 199(852) TR-X12	 2.2(1.3)/11.0 29.0 200(852L) TR-X15	 1.4/8.0 28.0 249(862) FO-X21	 1.6/8.0 28.0 249(862) FO-X22			
173-X014M	173-X016M	173-X018M	199-X016M	200-X022M	249-X014M	249-X016M			

Super Coarse

 1.6 25.0 001(801) BR-X40S	 1.8 25.0 001(801) BR-X31S			 1.4 25.0 001(801) BR-X41C	 1.6 25.0 001(801) BR-X40C	 1.8 25.0 001(801) BR-X31C
001-X016S	001-X018S			001-X014C	001-X016C	001-X018C

RA,CA 204 22.0mm

 1.2 22.0 001(801) BR-R46	 1.4 22.0 001(801) BR-R41	 1.6 22.0 001(801) BR-R40	 1.8 22.0 001(801) BR-R31	 2.3 22.0 001(801) BR-R26	 1.0(0.8)/1.0 22.0 010(805) SI-R45	 1.2(0.9)/1.2 22.0 010(805) SI-R46	 1.4(1.0)/1.6 22.0 010(805) SI-R47	 1.6(1.1)/1.8 22.0 010(805) SI-R48	 1.0/4.0 22.0 109(836) SF-R41
001-R012M	001-R014M	001-R016M	001-R018M	001-R023M	010-R010M	010-R012M	010-R014M	010-R016M	109-R010M
 1.2/4.0 22.0 109(836) SF-R31	 1.4/4.0 22.0 109(836) SF-R42	 1.2/8.0 24.0 111(837L) SF-R11	 1.4/8.0 24.0 111(837L) SF-R12	 1.6/8.0 24.0 111(837L) SF-R13	 1.2/8.0 24.0 141(881) SR-R11	 1.4/8.0 24.0 141(881) SR-R12	 1.4/8.0 24.0 165(858) TC-R21	 1.6/10.0 26.0 166(859) TC-R11	 1.0(0.6)/3.5 22.0 169(845) TF-R41
109-R012M	109-R014M	111-R012M	111-R014M	111-R016M	141-R012M	141-R014M	165-R014M	166-R016M	169-R010M
 1.2(0.7)/4.0 22.0 170(845) TF-R42	 1.8(1.0)/8.0 25.0 198(850) TR-R23	 1.6(0.8)/10.0 26.5 199(852) TR-R12	 1.8/3.6 22.0 257(368,830) FO-R32						
170-R012M	198-R018M	199-R016M	257-R018M						

Fine

 1.6 22.0 001(801) BR-R40F	 1.8 22.0 001(801) BR-R31F	 2.3 22.0 001(801) BR-R26F	 1.4/8.0 24.0 165(858) TC-R21F	 1.8/3.6 22.0 257(368,830) FO-R32F				 1.6 22.0 001(801) BR-R40EF	 1.8 22.0 001(801) BR-R31EF	 2.3 22.0 001(801) BR-R26EF
001-R016F	001-R018F	001-R023F	165-R014F	257-R018F				001-R016EF	001-R018EF	001-R023EF


Extra Fine



Tungsten Carbide Bur

Advantages

- 1 . Strictly comply with ISO Standard
- 2 . Sharper and More Efficient
- 3 . Durable ond Voluable
- 4 . FG, FG Short, FG Long, RA, RA Long etc



Tungsten Carbide Bur

Ball, Plain Cut	FG 1/4	FG 1/2	FG 1	FG 2	FG 3	FG 4	FG 5	FG 6	FG 7
FG 8	FG 9	Inverted Cone	FG 34	FG 35	FG 36	FG 37	FG 38	FG 39	Cylinder, Plain Cut
FG 56	FG 57	FG 58	FG 59	FG 60	Flat End Taper, Plain Cut	FG 168	FG 169	FG 170	FG 171
FG 172	FG 169L	FG 170L	FG 171L	FG 172L	Pear Plain Cut	FG 245	FG 246	FG 329	FG 330
FG 331	FG 332	FG 331L	FG 332L	X-Mas Tree, 3 slided	Carving of tooth pattern (Morphology)	FG 468	Cylinder, Cross Cut	FG 555	FG 556
FG 557	FG 558	FG 559	FG 560	FG 556L	FG 557L	FG 558L	FG 559L	Flat End Taper, Cross Cut	FG 698

FG 699	FG 700	FG 701	FG 702	FG 699L	FG 700L	FG 701L	FG 702L	End Cutting	FG 957
FG 958	Round Cylinder, Plain Cur	FG 1155	FG 1156	FG 1157	FG 1158	FG 1159	FG 1160	Round Taper, Plain Cut	
FG 1169	FG 1170	FG 1171	FG 1172	Round Cylinder, Cross Cut	FG 1556	FG 1557	FG 1558	FG 1559	
Round Taper, Cross Cut	FG 1701	FG 1702	Pear, Cross Cut	FG 1931	FG 1932	Round Cyliner, Cross Cut	FG 1957		
FG 1958	Round Taper, Cross Cut	FG 2057	FG 2058	Scaling to remove Tartar	ET3	ET5			

FG Short 313 16.5mm

FG S1	FG S2	FG S3	FG S4	FG S5	FG S245	FG S330	FG S331	FG S556	FG S557	FG S1157






FG Surgical 316 25mm/1 inch

 0.8 25.0 001	 1.0 25.0 001	 1.2 25.0 001	 1.4 25.0 001	 1.6 25.0 001	 1.8 25.0 001	 2.3 25.0 001	 1.0/4.0 25.0 107	 1.2/4.0 25.0 107	 1.0/4.0 25.0 168
FG X1	FG X2	FG X3	FG X4	FG X5	FG X6	FG X8	FG X57	FG X58	FG X170
 1.2/4.0 25.0 168	 1.0/4.0 25.0 107	 1.2/4.0 25.0 107	 1.4/4.0 25.0 107	 1.0/4.0 25.0 168	 1.2/4.0 25.0 168	 1.6/4.0 25.0 168	 1.0/4.0 25.0 137	 1.2/4.0 25.0 137	
FG X171	FG X557	FG X558	FG X559	FG X700	FG X701	FG X702	FG X1557	FG X1558	

RA,CA 204 22.0mm

 0.6 001	 0.8 001	 1.0 001	 1.2 001	 1.4 001	 1.6 001	 1.8 001	 2.1 001	 2.3 001	 0.8 010
RA 1/2	RA 1	RA 2	RA 3	RA 4	RA 5	RA 6	RA 7	RA 8	RA 34
 1.0 010	 1.2 010	 1.4 010	 1.6 010	 1.8 010	 2.1 010	 2.3 010	 0.9/2.5 233	 0.9/4.0 107	 1.0/4.0 107
RA 35	RA 36	RA 37	RA 38	RA 39	RA 40	RA 41	RA 245	RA 556	RA 557
 1.2/4.0 107	 1.4/4.0 107	 1.6/4.0 107	 1.8/4.0 107	 2.1/4.0 107	 2.3/4.0 107	 0.9/4.0 168	 1.0/4.0 168	 1.2/4.0 168	 1.6/4.0 168
RA 558	RA 559	RA 560	RA 561	RA 562	RA 563	RA 699	RA 700	RA 701	RA 702
 2.1/4.0 168	 2.3/4.0 168								
RA 703	RA 704								

RA,CA Long205 26mm

 1.0 001	 1.2 001	 1.4 001	 1.6 001	 1.8 001	 2.1 001	 2.3 001	 1.0/4.0 107	 1.2/4.0 168	
RA X2	RA X3	RA X4	RA X5	RA X6	RA X7	RA X8	RA X557	RA X701	
 1.6(0.9)/11 23.0 199	 1.6(0.9)/11 25.0 199	 1.6(0.9)/11 28.0 199		 1.6(1.0)/10.0 25.0 408		 1.6(1.0)/10.0 25.0 408	 1.2(0.7)/6.0 28.0 415	 1.2(0.7)/6.0 34.0 415	
Round End Taper	151	X151	XL151	Bone Cutter	161	Bone Cutter	162	H253E	H254E
 1.6(1.2)/8.0 25.0 172	 1.6(1.2)/8.0 28.0 172	 1.6(1.2)/8.0 32.0 172	 1.6/8.0 25.0 111	 1.6/8.0 28.0 111		 1.6(0.9)/9.0 21.0 210	 1.6(0.9)/9.0 23.0 210	 1.6(0.9)/9.0 25.0 210	
Flat End Taper,Cylinder	EC25	EC28	EC32	ES25	ES28	Round End Taper	152	X152	XL152

Finishing Carbide Burs

 1.4 19.0 001	 1.6 19.0 001	 1.8 19.0 001		 1.4/9.0 23.0 184	 1.6/9.0 23.0 184		 1.0/4.2 19.0 159	 1.4/6.0 21.0 164	
Round	FG 7004	FG 7005	FG 7006	Flat End Taper	FG 7204	FG 7205	Pointed Taper	FG 7611	FG 7612
 1.4/3.3 19.0 277	 1.8/3.5 19.0 277		 0.9/3.5 19.0 168	 1.0/5.5 21.0 170	 1.4/9.0 23.0 172		 1.6(1.0)/8.0 22.0 198		
Egg	FG 7404	FG 7406	Flat End Taper	FG 7114	FG 7214	FG 7714	Round End Taper	FG 7675	
		 0.9/3.6 19.0 496	 1.0/3.6 19.0 496	 1.2/3.6 19.0 496	 1.2/9.0 23.0 250				
8-12 Blades	Flame	FG 7901	FG 7902	FG 7903	FG 7903L				

IPR Set




BUR004

Veneer Set



BUR005
EX-31、FO-30EF、FO-32C、FO-41F、SF-11
TC-21F、TR-13C、TR-13F、TR-62F、VP-20


Crown Preparation Set



BUR006
CR-11F、FO-30F、FO-32C、FO-54C、SF-11、SF-12C、TR-13C、TR-62C、TR-62F

10 Holes Dental Burs

10 Holes Dental Burs Block Holder Aluminum



BUR001

24 Holes Dental Burs

24 Holes Dental Burs Block Holder Aluminum



BUR002

Lab Bur

C1517K



C1517K
Packing:1 pc/per pack

H0615PS



H0615PS
Packing:1 pc/per pack

F0615M



F0615M
Packing:1 pc/per pack

SD-06S



SD-06S
Packing:1 pc/per pack

L358S



L358S
Packing:1 pc/per pack

L0208XC



L0208XC
Packing:1 pc/per pack

C0208S



C0208S
Packing:1 pc/per pack

D040S



D040S
Packing:1 pc/per pack

AMD0918D



AMD0918D
Packing:1 pc/per pack

C0613



C0613
Packing:1 pc/per pack

D-4



D-4
Packing:1 pc/per pack

SJ-08S



SJ-08S
Packing:1 pc/per pack

Double Stainless Steel Instrument Trolley



808-004
Shelf Capacity: 80kg
Overall Size (WxDxH): 800*500*900mm
Unit Weight: 14kg

5 Tiers Medical CliniCart Instrument Trolley



Size : 480 * 505 * 820
Stainless steel
5 tier drawers with 2 tiers of instruments division compartments
Safety locks
Slide in-out easy working station
Balanced wheels with locks
Soft closing drawers
Smooth welding
Easy mobility



808-003

Medical Motor CliniCart With Build-In Power Outlet

Size : 475 * 450 * 825
Stainless steel
Build-in power outlet for surgical motors
2 tier drawers
Balanced wheels with locks
Soft closing drawers
Smooth welding
Easy mobility



808-002

DTE Implant motor



DTE-1
Speed Range: 300- 40,000rpm
Torque Range (20:1): 5N.cm- 80N.cm
The Motor Weight: 140g
The Largest Water: 110ml/min
Power Supply: AC220V 50/60Hz
Volume: Mute
1 Year warrenty



Outstanding quality makes it your solid choice

Implant-X adopts mute motor imported from Swiss, Advanced FOC control ng technique blends low speed performance with highspeed performance, In order to verify the clinical features of flexible response and stable safety, Implant-XR&D personnels go into the operating room, getting them involved in clinical applications, Under the instruc ion and assistance of experts from home and abroad, they finally complete the clinical optimization of product, DTE spares no efforts to provide you high quality products.

Intuitive userinterface

Large color screen enables clear display, Visual images display every procedure in implant surgery, Preset parameters eliminate setting process.

Specification

Rotation Speed Range	300-40.000rpm
Torque Range/Gear 20: 1	5N.cm-80N.cm
Max Pump Output	110ml/min
Main Unit Dim nen sions	276mm*267mm*110mm
Motor Dimensions	φ23.4mm*L91mm
Motor Weight	140g
PowerSupply	AC220V50/60Hz



Strong and stable contra-angle

LED contra-angle ensures visibility, helping to achieve treatment of high accuracy and eff ciency. High precision contra-angle with tiny oscillation enables stable torque output under highspeed and low speed. Imported advanced bearings guarantee performance and durability of contra-angle.



Color touchscreen with glass surface

The colour touchscreen with its glass surface makes operation of the control unit simple and intuitive, The glass surface is easy to clean Wiping disinfection ensures safety in treatment area



Mute peristaltic pump

Mute peristaltic pump water supply system is easy to connect Be compatible with disposable oral water pipe to ensure safety of treatment.



Flexible multi-function foot pedal

Muli-function foot pedal can realize aseptic operation and liberate two han

Baxter Sodium Chloride 0.9% for Irrigation



OS7124
Specification: 1L
Packing: Bottle



OS1323
Specification: 500ml
Packing: Bag

Kiato Surgical Blades

Individually packed in easy peel color coded foil & boxes
Sterilised by Gamma Radiation to a minimum dosage of 2.5 MRads



OSK12
Packing: Box of 100 blades



OSK15
Packing: Box of 100 blades




Cardinal Health Protexis Latex Classic Powder Free Surgical Gloves



SKU	Specification	Packing
OSN55X	Cream/ Light Brown – Size 5.5	50 pcs/box
OSN60X	Cream/ Light Brown – Size 6.0	50 pcs/box
OSN70X	Cream/ Light Brown – Size 7.0	50 pcs/box
OSN75X	Cream/ Light Brown – Size 7.5	50 pcs/box
OSN80X	Cream/ Light Brown – Size 8.0	50 pcs/box
OSN85X	Cream/ Light Brown – Size 8.5	50 pcs/box
OSN90X	Cream/ Light Brown – Size 9.0	50 pcs/box

Cardinal Health Protexis® Micro PI Latex-Free



SKU	Specification	Packing
OSPM55	Size – 5.5	50 pcs / box
OSPM60	Size – 6.0	50 pcs / box
OSPM65	Size – 6.5	50 pcs / box
OSPM70	Size – 7.0	50 pcs / box
OSPM75	Size – 7.5	50 pcs / box
OSPM80	Size – 8.0	50 pcs / box
OSPM85	Size – 8.5	50 pcs / box
OSPM90	Size – 9.0	50 pcs / box

Cardinal Health Protexis® Neoprene Essential



SKU	Specification	Packing
OSDS55	Size – 5.5	50 pcs / box
OSDS60	Size – 6.0	50 pcs / box
OSDS65	Size – 6.5	50 pcs / box
OSDS70	Size – 7.0	50 pcs / box
OSDS75	Size – 7.5	50 pcs / box
OSDS80	Size – 8.0	50 pcs / box
OSDS85	Size – 8.5	50 pcs / box
OSDS90	Size – 9.0	50 pcs / box

Cardinal Health Protexis® PI Latex-Free



SKU	Specification	Packing
OSPT55X	Size – 5.5	50 pcs / box
OSPT60X	Size – 6.0	50 pcs / box
OSPT65X	Size – 6.5	50 pcs / box
OSPT70X	Size – 7.0	50 pcs / box
OSPT75X	Size – 7.5	50 pcs / box
OSPT80X	Size – 8.0	50 pcs / box
OSPT85X	Size – 8.5	50 pcs / box
OSPT90X	Size – 9.0	50 pcs / box

Non - Absorbable suture

LOOK Silk Suture; Black Braided;18in/45cm;Reverse Cutting;19mm;3/8 Circle



SKU	Specification	Packing
OS784B	5-0; 19mm	12 pcs/Box
OS781B	4-0; 19mm	12 pcs/Box
OS754B	3-0; 19mm	12 pcs/Box

LOOK Nylon; Black Monofilament; Reverse Cutting



SKU	Specification	Packing
OS928B	3-0; 19mm	12 pcs/Box
OS922B	4-0; 19mm	12 pcs/Box
OS920B	4-0; 24mm	12 pcs/Box
OS918B	5-0; 19mm	12 pcs/Box
OS919B	5-0; 16mm	12 pcs/Box

LOOK PTFE Suture;19mm;3/8 Circle



SKU	Specification	Packing
OS816B	C-6 undyed monofilament; 3-0; 18"/45cm; reverse cutting; 18mm; 3/8 circle	12 pcs/Box
OS821B	C-3 undyed monofilament; 6-0; 18"/45cm; precision reverse cutting; 13mm; 3/8 circle	12 pcs/Box
OS822B	PC-31 undyed monofilament; 4-0; 18"/45cm; precision reverse cutting; 13mm; 3/8 circle	12 pcs/Box
OS823B	C-3 undyed monofilament; 5-0; 18"/45cm; precision reverse cutting; 13mm; 3/8 circle	12 pcs/Box

Absorbable suture

PolySyn FA™ BRAIDED FAST ABSORBABLE SUTURE; Undyed; DSM 19; 3/8 Circle Precision Reverse Cutting 19mm



SKU	Specification	Packing
OSF426	4-0; undyed, 27"	12 pcs/Box
OSF496	4-0; undyed, 18"	12 pcs/Box
OSF497	3-0; undyed, 18"	12 pcs/Box

PolySyn™ BRAIDED ABSORBABLE SUTURE; 19mm



SKU	Specification	Packing
OS421B	5-0; 27in/70cm; reverse cutting 19mm; 3/8 Circle	12 pcs/Box
OS422B	4-0; 27in/70cm; reverse cutting 19mm; 3/8 Circle	12 pcs/Box
OS423B	3-0; 27in/70cm; reverse cutting 19mm; 3/8 Circle	12 pcs/Box
OS422N	4-0; undyed, 27", DSM19, 3/8 circle reverse cutting 19mm	12 pcs/Box
OS426N	4-0; undyed, 27", DSM19, 3/8 circle precision reverse cutting 19mm	12 pcs/Box
OS423N	3-0; undyed, 27", DSM19, 3/8 circle reverse cutting 19mm	12 pcs/Box

Monoderm Monofilament Absorbable Suture; DSM19, 3/8 Circle Precision Reverse Cutting 19mm



SKU	Specification	Packing
OS495N	5-0; undyed, 18"	12 pcs/Box
OSF496	4-0; undyed, 18"	12 pcs/Box
OS497N	3-0; undyed, 18"	12 pcs/Box

Tooth Forceps For Adults

Multihole handle



652A
54|45
For upper bicuspid



650A
321|123
For upper Incisors



651A
321|123
For lower Incisors



653A
54|45
For lower bicuspid



654A
76|
For upper molars, right



655A
|67
For upper molars, left



657A
876|678
For lower molars

Tooth Forceps For Adults

Concave handle



665C
For upper and lower roots



665C-1
For upper and lower residual roots



663C-1
For upper roots



664C-1
For lower roots



659C-1
For the horizontal
embedded eighth molars

Tooth Forceps For Adults

Multihole handle



658A
8|8
For upper third molars



660A
76|
For upper right molars or residual crown



661A
|76
For upper left molars or residual crown



662A
76|67
For lower molars or residual crown



663A
▲▲
For upper and lower roots



665A
▲▲
For upper and lower roots



665E-2
▲▲
For upper and lower residual roots or fragments



665E-3
▲▲
For upper and lower residual roots or fragments

Tooth Forceps For Children



681A
321|123
For upper incisors



682A
321|123
For lower incisors



683A
54|45
54|45
For upper and lower deciduous molars



684A
▲▲
For upper and lower roots



685C
321|123
321|123
For deciduous laterals



686C
54|45
For upper deciduous molars



687C
54|45
For lower deciduous molars

Tooth Forceps For Children



681A
321|123
For upper incisors



682A
321|123
For lower incisors



683A
54|45
For upper and lower deciduous molars



684A
54|45
For upper and lower roots



685C
321|123
For deciduous laterals



686C
54|45
For upper deciduous molars



687C
54|45
For lower deciduous molars

Tooth Forceps



656B
anterior forceps



655A
Upper Left Molar Forceps



659B
Upper Incisor Root Forceps



673A
Lower molar wisdom forceps



662B
Upper Wisdom Forceps



657B
Upper Molar Forceps




654A
Upper Right Molar Forceps




SET-B7

Delicate Root Elevators



Delicate Root Elevators



Hand Instrument

Delicate Root Elevator Cryer



710-121C

710-122C

Delicate Root Elevator Root Picker



710-125C 5#
2mm
for root apex left-facing curved

710-124C 4#
2mm
for root apex straight

710-126C 6#
2mm
for root apex right-facing curved

Hand Instrument

Delicate Root Elevators

Pen-Like Handle



1#
710-401A 2.5mm
straight blade

2#
710-402A 4.5mm
straight blade

3#
710-403A 2.5mm
curved blade

4#
710-404A 4.5mm
curved blade

5#
710-405A 3.8mm
triangle blade

6#
710-406A 1.5/2.5mm
double blade

Trephine Burs






818-101A 1#A Φ3.5mm
818-102A 2#A Φ4.0mm
818-103A 3#A Φ4.5mm
818-104A 4#A Φ5.0mm
818-105A 5#A Φ5.5mm
818-106A 6#A Φ6.0mm
818-107A 7#A Φ6.5mm
818-108A 8#A Φ7.0mm

818-201A 1#A Φ3.0mm
818-202A 2#A Φ4.0mm
818-203A 3#A Φ5.0mm

818-206A 6#A Φ3.5mm
818-207A 7#A Φ4.0mm
818-208A 8#A Φ4.5mm
818-209A 9#A Φ5.0mm

SET- 821 Implant Universal Surgical Kit



Root elevator – Vertical lifter type I



710-611 1#
710-614 4# Niti end



710-612 2#
710-615 5# Niti end



710-613 3#
710-616 6# Niti end

Implant Screws



821-101L 1# (1.0)
821-103L 3# (1.27)
821-105L 5# (NOB)
821-107L 7# (0.9)
821-109L 9# (1.25)



821-102L 2# (1.2)
821-104L 4# (1.4)
821-106L 6# (ITI)
821-108L 8# (1.16)
821-110L 10# (TM)



821-101S 1# (1.0)
821-103S 3# (1.27)
821-105S 5# (NOB)
821-107S 7# (0.9)
821-109S 9# (1.25)



821-102S 2# (1.2)
821-104S 4# (1.4)
821-106S 6# (ITI)
821-108S 8# (1.16)
821-110S 10# (TM)

Universal Implant Driver Kit



820-101 1#



821-001

Tissue Scissors



705-112 14cm
705-132 16cm
curved/blunt,serrated



705-111 14cm
705-131 16cm
curved/sharp serrated



705-119 14cm
705-139 16cm
curved serrated



705-114 14cm
705-134 16cm
straight/sharp
705-116 14cm
705-136 16cm
straight/sharp serrated



705-115 14cm
705-135 16cm
straight/blunt
705-117 14cm
705-137 16cm
straight/blunt serrated



705-114B 12cm
straight/sharp serrated



705-111B 12cm
curved serrated

Needle Holders With TC



704-121 14cm
704-131 16cm
straight,narrow
704-122 14cm
704-132 16cm
straight,slender
704-125 14cm
704-135 16cm
straight,wide
704-134 16cm
curved,slender
704-322 14cm
704-352 17cm
straight,slender
curved handle

Suture Scissors



744-102 2#
744-201 1#

Micro Needle Holders



791-223 16cm
type, angled 90°
791-111 14cm
791-121 16cm
791-131 18cm
tissue, straight
791-112 14cm
791-122 16cm
791-132 18cm
tissue, curved

Micro Scissors




792-121 16cm
792-131 18cm
straight

792-122 16cm
792-132 18cm
curved

792-211 14cm
792-221 16cm
792-231 18cm
straight

792-212 14cm
792-222 16cm
792-232 18cm
curved

Micro Tweezers



790-111 14cm
790-121 16cm
790-131 18cm
straight

790-122 16cm
straight, single tooth

790-123 16cm
790-133 18cm
curved

790-124 16cm
curved, single tooth

Mosquito Forceps



739-112 12cm
739-122 14cm
curved



739-111 12cm
739-121 14cm
straight

Towel forceps



731-112 12cm
731-122 14cm
curved

Towel tweezer



741-101 1#

Scalpel Handles



1# 703-001 with elevator
2# 703-002
3# 703-003
4# 703-004
5# 703-005 movable
7# 703-007 automatic disassembly

Hand Instrument

Fine Tweezers



702-111 0.9mm
702-112 1.2mm
12cm
single tooth



702-211 0.9mm
702-212 1.2mm
12cm
single tooth



702-123 1.5mm
18cm
single tooth



702-124
18cm
multi-teeth



702-224
18cm
multi-teeth

Tweezers



702-115 12cm
702-135 15cm
flat teeth



702-136 15cm
scatheless



702-325 18cm
flat teeth

Suture Tweezers



805-101

Hand Instrument

Cheek Retractors



518-501 14X17mm 1#
518-502 20X22mm 2#
518-503 26X30mm 3#



518-511
4#



518-701
movable tip



518-301



518-302



518-303



518-204
double tips

Bone Calipers



708-001 17cm
inward



708-002 18cm
straight



708-003 18cm
curved



708-004 9cm
straight



708-005 9cm
curved



708-006 15cm
bone meters

Dental implant Guidance Appliance



816-101 1#




816-102 2#

Gag Gauges



708-201 1#



708-202 2#



Bone Mills



727-101
plier-like



727-201 1#
simple



727-202 2#
simple

Bone Placing Container



814-101 1# Φ38x28
814-102 2# Φ25x18



814-201 1# Φ38x28
with cover

Hand Instrument

Bone Carriers



735-101 1# Φ 2.5mm
735-102 2# Φ 3.5mm
735-103 3# Φ 4.5mm
735-104 4# Φ 6.0mm



735-201 1#

Bone Scraper



815-101
straight


815-102
urved



807-101 Bone Graft Scoop



807-201 1# Bone Lifter



807-202 2# Bone Lifter

Hand Instrument

Suction Tubes



734-711 Φ 1.5mm
Titanium tip



734-701 Φ 1.5mm
734-702 Φ 3.0mm
Titanium tip



734-301



734-201
with filter



734-601 light
734-101 heavy



734-501 heavy

Hand Instrument

Mallets



Mouth Mirror Anti-heat Plain



808-102A
Packing: 12pcs/box

Mouth Mirror Handle

Standard dental mirror handle for use with dental mirror.



808-101A
Packing: 1pcs/bag

Hand Instrument

Dental Retractors



819-100A




819-101A 1# 27.5mm
819-102A 2# 15.5mm
819-103A 3# 6.5mm
819-104A 4# 6.5mm
819-105A 5# 6.5mm




819-201B
1#


Periosteal Elevators



Bone Files




Jowl Chisels



Jowl Chisels


round handle, single slanting cutting edge (with scale)




round handle, double slanting cutting edge (with scale)



Enamel Chisel




Surgical Curettes




Bone Spreading System

803-100 8X1 (Set)




Sinus Lifts



Hand Instrument

Periodontal Files



Surgical Curette



Hand Instrument

SET-818 Bone Collection Kit




Hand Instrument

SET-820 Ridge Splitting Kit



Hand Instrument



Gingival Separators



Hand Instrument

Plastic Plugger


807-310A





Hand Instrument


Plastic Plugger

807-310F






Rubber Dam Clamp




Product number	Specification	purpose
740-111	Q1#	For centrals
740-112	Q2#	For anteriors
740-121	S1#	For small bicuspid and pedo central
740-122	S2#	or Lower bicuspid with flat beaks
740-123	S3#	For bicuspid, but with larger beaks
740-124	S4#	For upper & lower bicuspid with festooned beaks
740-125	S5#	For upper & lower bicuspid with flat beaks
740-131	M1#	For upper molars
740-132	M2#	For molars but with festooned beaks
740-133	M3#	For Lower molars with flat beaks
740-134	M4#	For Lower molars with festooned beaks
740-135	M5#	For Large molars with larger beaks
740-141	C1#	For lower molars/wisdom teeth, serrated
740-142	C2#	For lower molars/wisdom teeth, serrated
740-143	C3#	For upper left and lower right molar, serrated
740-144	C4#	For upper right and lower left molar, serrated


Periodontal Chisels




Canal Plugger




Gingival Retraction Cord Plugger



Dental Probes



Periodontal Knives



Basic Implant Instrument Kit



SET-1




Micro Surgical Suture Kit



SET-2



General Surgical Suture Kit



SET-3



Ortho Mirrors



500-300

Ultra Supersonic Technology Powered Toothbrush

Are clinically proven to be effective against gum disease. Dental Experts from Academic, Hospitals and Industry strongly recommended these powered toothbrushes. Daily use has been proven to be significant reduce oral bacteria counts and improve gingival regeneration.

- Superb Cleaning Power
- Small Head And Soft Bristles
- Wireless Charging & Long Battery Life
- Slim Body For Easy Grip
- Waterproof IPX7
- Interchangeable Head- Has Facial Cleaning And Massager Function



IC038
Specification: Orange Yellow



IC039
Specification: Green

Tooth Model Soft



S001

Tooth Model Hard



S002